

STANDARD DEVICE

AE 1+2

Abrasion and friability testers

The abrasion and friability testing device consists of two test stations. If required, test drums can be equipped to test for wear due to rolling and falling or to test the abrasion of tablets, tablet cores, granules etc.

The tester meets all the latest requirements of the European & US Pharmacopoeia, GMP and GLP.

- · Quick drum assembly
- · User-friendly membrane keypad
- GMP-compliant stainless steel housing

MODEL	LxWxH(mm)	Weight (kg)
AE1	320 x 300 x 250	6.5
AE 2	320 x 430 x 250	7

MEASURED PARAMETERS



ABRASION



FRIABILITY

ADVANTAGES

- + Compact stainless steel housing
- + Quick drum assembly
- + Automatic emptying

OPTIONS

- · Friability test drum
- · Abrasion test drum





Abrasion and friability testing device AE 2 with 2 test drums